



NC807C

Ceiling Pull Unit



Overview

Can be used to generate a standard call on an 800 Series Conventional Call System.

Typically wired as a slave to an 800 Series Call Point or Reset Point so that its calls can be cancelled.

Supplied with 3m of red nylon pull cord and two red open-sided triangular bangles for ease of use by the infirm.

Red confidence LED and diffuser emits a bright red glow over a wide area to reassure the user that a call has been made.

Surface mountable.

Base fits 16mm conduit and BESA centres.

NCP-13

Anti-Bacterial Wipe Clean Pull Cord Accessory Pack



Overview

A BS8300 compliant anti-bacterial wipe-clean pull cord pack suitable for use with C-TEC's Quantec and 800 Series ceiling pull units.

Contains 3m of red vinyl encapsulated cord, two triangular bangles and a red fixing bullet.

Cord has integral properties that can help prevent the spread of infections such as MRSA, e.coli, c.diff, S Aureus and the Norovirus.

Ideal for use in hospitals, nursing and care homes, doctors' surgeries, clinics and any other area where infection control is imperative.

A simple and cost-effective means of updating existing ceiling pulls to comply with the stringent hygiene regulations of many health authorities and inspection bodies.

Easy to clean using a clean cloth and mild detergent.

Also compatible with many other manufacturers' ceiling pull units.

Replacement nylon pull cord packs also available (part no. NCP-12)

For more information on these products, please contact us:



T: 00353 47 86274 / 00353 47 86663

E: Sales@irishsafetysystems.ie